

Anti-DLAT antibody (1-290) (STJ28210)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IF/CC:1:20-1:50 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 69kDa Observed Mw: 69kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	1737
Gene Symbol	DLAT
UniProt ID	ODP2_HUMAN
Immunogen Region	1-290
Immunogen Sequence	MWRVCARRAQNVPWAGLEA RWTALQEVPGTPRVTSRSGP APARRNSVTTGYGGVRLCG WTPSSGATPRNRLQLLGS PGRYYSLPPHQKVPLPSLS PTMQAGTIARWEKKEGDKIN EGDLIAEVETDKATVGFESL EECYMAKILVAEGTRDVPIG AIICITVGKPEDIEAFKNYT LDSSAAPTQAAPAPTPAAT ASPPTPSAQAPGSSYPHMQ VLLPALSPMTMTMGTVQRWE
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human DLAT (NP_001922.2).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunofluorescence analysis of HeLa cells using DLAT antibody (STJ28210). Blue: DAPI for nuclear staining.



Western blot analysis of extracts of various cell lines, using DLAT antibody (STJ28210) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.