

Anti-ST3GAL4 antibody (27-267) (STJ28231)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:100-1:500 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 38kDa Observed Mw: 50-55kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6484
Gene Symbol	ST3GAL4
UniProt ID	SIA4C_HUMAN
Immunogen Region	27-267
Immunogen Sequence	REDSFYFPIPEKKEPCLQGE AESKASKLFGNYSRDQPIFL RLEDYFWVKTPSAYELPYGT KGSEDLLRLVLAITSSSIPK NIQLRCRRCVVVGNGHRLR NSSLGDAINKYDVVIRLNNA PVAGYEGDVGSKTTMLRFYP ESAHFDPKVENNPDTLLVLV AFKAMDFHWIETILSDKKRV RKGFWKQPPLIWDVNPQIR ILNPFMEIAADKLLSLPMQ QPRKIKQKPTTGLLAITLA
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 27-267 of human ST3GAL4 (NP_006269.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of Mouse large intestine, using ST3GAL4 antibody (STJ28231) at 1:400 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.