

Anti-TDG antibody (1-410) (STJ28323)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 46kDa Observed Mw: 50kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6996
Gene Symbol	TDG
UniProt ID	TDG_HUMAN
Immunogen Region	1-410
Immunogen Sequence	MEAENAGSYSLQQAQAFYTF PFQQLMAEAPNMAVVNEQQM PEEVPAPAPAEQEPVQEAPKG RKRKPRTTEPKQPVEPKPV ESKKSGKSAKSKEKQEKTD TFKVKRKVDRFNGVSEAELL TKTLPDILTFNLDIVIIGIN PGLMAAYKGHHYPGPGNHFH
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human TDG (NP_003202.3).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of extracts of Rat thymus, using TDG antibody (STJ28323) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.