

Anti-NFKBIB antibody (1-356) (STJ28344)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/IHC-P/ELISA
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	Lot specific
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	WB:1:500-1:2000 IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
	Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
	Isotype	IgG
	Molecular Weight	Protein Mw: 38kDa Observed Mw: 40-50kDa
	Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	4793
	Gene Symbol	NFKBIB
	UniProt ID	IKBB_HUMAN
	Immunogen Region	1-356
	Immunogen Sequence	MAGVACLGKAADADEWCDSG LGS LGPDAAAPGGPGLGAEL GPGLSWAPLVFGYVTEGDGT ALHLAVIHQHEPFLDFLLGF SAGTEYMDLQNDLQQTALHL AAILGETSTVEKLYAAGAGL CVAERRGHTALHLACRVGAH ACARALLQPRPRRPREAPDT YLAQGPDRTPDNTNHTPVLY PDSLDKEEEEESEEDWKLQL EAENYEGHTPLHVAVIHKDV EMVRLLRDAGADLDKPEPT
	Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-356 of human NFKBIB (NP_002494.2).
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
   		