

Anti-SMARCAD1 antibody (1-300) (STJ28413)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IF/CC:1:50-1:100 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 117kDa Observed Mw: Refer to figures
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	56916
Gene Symbol	SMARCAD1
UniProt ID	SMRCD_HUMAN
Immunogen Region	1-300
Immunogen Sequence	MNLFNLD RFRFEKRNKIEEA PEATQPSQPGPSSPISLSA EEENAEGEVSRANTPDSDIT EKTEDSSVPETPDNERKASI SYFKNQRGIQYIDLSSDSED VVSPNCSNTVQEKTFNKDTV IIVSEPSDEESQGLPTMAR RND DISELEDLSELEDL KDA KLQTLKELFPQRSDNDLLKL IESTSTMDGAIAAALLMFGD AGGGPRKRKLSSSSSEPYEED EFND DQSIKKTRLDHGEEES
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human SMARCAD1 (NP_001121901.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Immunofluorescence analysis of A549 cells using SMARCAD1 antibody (STJ28413). Blue: DAPI for nuclear staining.