

Anti-CLIC1 antibody (42-241) (STJ28446)

GENERAL INFORMATION

	Product Type	Primary antibodies
	Applications	WB/ELISA
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 27kDa Observed Mw: 27kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1192
Gene Symbol	CLIC1
UniProt ID	CLIC1_HUMAN
Immunogen Region	42-241
Immunogen Sequence	NVTTVDTKRRTETVQKLCPG GQLPFLLYGTEVHTDTNKE EFLEAVLCPPPRYPKLAALNP ESNTAGLDIFAKFSAYIKNS NPALNDNLEKGLLKALKVLD NYLTSPLPEEVDETSAEDEG VSQRKFLDGNELTLADCNLL PKLHIVQVVCCKKYRGFTIPE AFRGVHRYLSNAYAREEFAS TCPDDEEIELAYEQVAKALK
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 42-241 of human CLIC1 (NP_001279.2).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of various lysates using CLIC1 Rabbit pAb (STJ28446) at 1:1000 dilution.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution.
Lysates/proteins: 25 μ g per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Enhanced Kit
Exposure time: 30s.