

Anti-PTPN3 antibody (1-300) (STJ28483)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 104kDa Observed Mw: 103kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5774
Gene Symbol	PTPN3
UniProt ID	PTN3_HUMAN
Immunogen Region	1-300
Immunogen Sequence	MTSRLRALGGRINNIRTSEL PKEKTRSEVICSIHFLDGVV QTFKVTQKDTGQVLLDMVHN HLGVTEKEYFGLQHDDSDVD SPRWLEASKAIRKQLKGGFP CTLHFRVRFFIPDPNTLQQE QTRHLYFLQLKMDICEGRLT CPLNSAVVLASYAVQSHFGD YNSSIHPGYLSDSHFIPDQ NEDFLTKEVSLHEQHSGLKQ SEAESCYINIARTLDFYGV E LHSGRDLHNLDLMIGIASA
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human PTPN3 (NP_001138840.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of extracts of various cell lines, using PTPN3 antibody (STJ28483) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 10s.