

Anti-CLEC3B antibody (22-202) (STJ28501)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF/ICC:1:10-1:100 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 21kDa Observed Mw: Refer to figures
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7123
Gene Symbol	CLEC3B
UniProt ID	TETN_HUMAN
Immunogen Region	22-202
Immunogen Sequence	EPPTQKPKKIVNAKKDVVNT KMFEELKSRLDTLAQEVALL KEQQALQTVCLKGTKVHMKC FLAFTQTKTFHEASEDCISR GGTLGTPQTGSENDALYEYL RQSVGNEAEIWLGLNDMAAE GTWVDMTGARIAYKNWETEI TAQPDGGKTENCAVLGAAN GKWFDKRCRDQLPYICQFGI V
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 22-202 of human CLEC3B/Tetranectin (NP_003269.2).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Immunofluorescence analysis of HeLa cells using CLEC3B/Tetranectin antibody (STJ28501). Blue: DAPI for nuclear staining.