

Anti-NGB antibody (1-151) (STJ28560)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA: Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 17kDa Observed Mw: 17kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	58157
Gene Symbol	NGB
UniProt ID	NGB_HUMAN
Immunogen Region	1-151
Immunogen Sequence	MERPEPELIRQSWRAVSRSP LEHGTVLFARLFALPDLLP LFQYNCRQFSSPEDCLSSPE FLDHIRKVMLVIDAAVTNVE DLSSLEEYLA SLGRKHRAVG VKLSSFSTVGESLLYMLEKC LGPAFTPATRAAWSQLYGAV VQAMSRGWDGE
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-151 of human NGB (NP_067080.1).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of lysates from Mouse brain, using NGB Rabbit polyclonal antibody (STJ28560) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.