

Anti-FRZB antibody (166-325) (STJ28675)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Rat/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:2000-1:7000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 36kDa Observed Mw: 36kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2487
Gene Symbol	FRZB
UniProt ID	SFRP3_HUMAN
Immunogen Region	166-325
Immunogen Sequence	SNGNCRGASSERCKCKPIRA TQKTYFRNNYNYVIRAKVKE IKTCHDVTAVVEVKEILKS SLVNIIPRDTVNLYTSSGCLC PPLNVNEEYIIMGYEDEERS RLLLVESGIAEKWKDRLGKK VKRWDMKLRHLGLSKSDSSN SDSTQSQKSGRNSNPRQARN
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 166-325 of human FRZB (NP_001454.2).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of lysates from Rat spleen, using FRZB Rabbit polyclonal antibody (STJ28675) at 1:6000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.