

## Anti-TDRKH antibody (36-355) (STJ28850)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/ELISA           |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                        |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                                                             |
| Concentration       | Lot specific                                                                                                                                           |
| Conjugation         | Unconjugated                                                                                                                                           |
| Purification        | Affinity purification                                                                                                                                  |
| Dilution Range      | WB:1:500-1:2000<br>ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements. |
| Formulation         | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.                                                                                                     |
| Isotype             | IgG                                                                                                                                                    |
| Molecular Weight    | Protein Mw: 62kDa<br>Observed Mw: 72kDa                                                                                                                |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                                         |

### TARGET INFORMATION

|                    |                                                                                                                                                                                                                                                                     |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">11022</a>                                                                                                                                                                                                                                               |
| Gene Symbol        | <a href="#">TDRKH</a>                                                                                                                                                                                                                                               |
| UniProt ID         | <a href="#">TDRKH_HUMAN</a>                                                                                                                                                                                                                                         |
| Immunogen Region   | 36-355                                                                                                                                                                                                                                                              |
| Immunogen Sequence | RYRESREERLTFVGEDDIEI EMRVPOEAVKLIIGRQGANI KQLRKQTGARIDVDTEDVGD<br>ERVLLISGFVPVQVCKAKAAI HQILTENTPVSEQLSVPQRS VGRIIGRGGETIRSICKASG<br>AKITCDKESEGTLLLSRLIK ISGTQKEVAAAKHLILEKVS EDEELRKRIHSAETRVPRK<br>QPISVRREDMTEPGGAGEPA LWKNTSSSMEPTAPLVTPPP KGGGDMAVVVSKEGSWEKP |
| Specificity        | Recombinant fusion protein containing a sequence corresponding to amino acids 36-355 of human TDRKH (NP_001077432.1).                                                                                                                                               |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|



Western blot analysis of extracts of various cell lines, using TDRKH antibody (STJ28850) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.