

Anti-LGALS3BP antibody (306-585) (STJ29009)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IF/CC:1:50-1:100 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 65kDa Observed Mw: 75kDa-110kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	3959
Gene Symbol	LGALS3BP
UniProt ID	LG3BP_HUMAN
Immunogen Region	306-585
Immunogen Sequence	LLPRSDLAVPSELALLKAV DTWSWGERASHEEVEGLVEK IRFPMMLPEELFELQFNLSL YWSHEALFQKKTQLQALEFHT VPFOLLARYKGLNLTEDTYK PRIYTSPTWSAFVTDSSWSA RKSQLVYQSRRGPLVKYSSD YFQAPSDYRYYPYQSFQTPQ HPSFLFQDKRVSWSLVYLP IQSCWNYGFSCSSDELPLVGLG LTKSGGSDRTIAYENKALML CEGLFVADVTDFEGWKAAI
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 306-585 of human LGALS3BP (NP_005558.1).

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Immunofluorescence analysis of MCF7 cells using LGALS3BP antibody (STJ29009). Blue: DAPI for nuclear staining.



Western blot analysis of extracts of various cell lines, using LGALS3BP antibody (STJ29009) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.