

## Anti-RAD51B antibody (1-384) (STJ29042)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Mouse/Rat

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 42kDa Observed Mw: 42kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">5890</a>
Gene Symbol	<a href="#">RAD51B</a>
UniProt ID	<a href="#">RA51B_HUMAN</a>
Immunogen Region	1-384
Immunogen Sequence	MGSKLKRVLGSLQELCDRLS RHQILTCQDFLCLSPLELMK VTGLSYRQVHELLCMVSRAC APKMQTAYGIKAQRSADFSP AFLSTTLSALDEALHGGVAC GSLTEITGPPGCGKTQFCIM MSILATLPTNMGGLEGAVVY IDTESAFSAERLVEIAESRF PRYFNTEEKLLLTSSKVHLY RELTCDVQLQRIESLEEEII SKGIKLVILDSVASVVRKEF DAQLQGNLKERNKFLAREA
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-384 of human RAD51B (NP_598193.2).

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Western blot analysis of various lysates using RAD51B Rabbit polyclonal antibody (STJ29042) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.