

Anti-CCNE2 antibody (255-404) (STJ29112)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IF/CC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 47kDa Observed Mw: 47kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	9134
Gene Symbol	CCNE2
UniProt ID	CCNE2_HUMAN
Immunogen Region	255-404
Immunogen Sequence	ALKDAPKVLPPQYSQETFIQ IAQLDLCILAIIDSLEFQYR ILTAAALCHFTSIEVVKKAS GLEWDSISECDWMVPFVNV VKSTSPVKLTKTFKKIPMEDR HNIQTHTNYLAMLEEVNYIN TFRKGGQLSPVCNGGIMTPP KSTEKPPGKH
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 255-404 of human Cyclin E2 (NP_477097.1).

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Immunofluorescence analysis of MCF7 cells using Cyclin E2 antibody (STJ29112).



Western blot analysis of extracts of Jurkat cells, using Cyclin E2 antibody (STJ29112). Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST.