



Anti-ICOSLG antibody (19-256) (STJ29160)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES	
Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:100-1:500 IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 33kDa Observed Mw: 60kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	23308
Gene Symbol	ICOSLG
UniProt ID	ICOSL_HUMAN
Immunogen Region	19-256
Immunogen Sequence	DTQKEVVRAMVGSDVELSCA CPEGSRFDLNDVYVYWTSE SKTVVTYHIPQNSSLENVDS RYRNRALMSPAGMLRGDFSL RLFNVTPQDEQKFHCLVLSQ SLGFQEVLSVEVTLHVAANF SVPVVSAPHSPSQDELTFCT TSINGYPRPNVYWKNTDNS LLDQALQNDTVFLNMRGLYD VSVLRIARTPSVNIGCCIE NVLLQQNLTVGSQTGNDIGE RDKITENPVSTGEKNAAT
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 19-256 of human ICOSL (NP_056074.1).
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunohistochemistry analysis of paraffin-embedded human kidney using ICOSL Antibody (STJ29160) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded human skin carcinoma using ICOSL Antibody (STJ29160) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.