

Anti-PCGF3 antibody (1-242) (STJ29595)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 28kDa Observed Mw: 28kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10336
Gene Symbol	PCGF3
UniProt ID	PCGF3_HUMAN
Immunogen Region	1-242
Immunogen Sequence	MLTRKIKLWDINAHITCRLC SGYLIDATTVTECLHTFCRS CLVKYLEENNTCPTCRIVIH QSHPLQYIGHDRMTMODIVYK LVPGLQEAEMRKQREFYHKL GMEVPGDIKGETCSAKQHLD SHRNGETKADDSSNKEAAEE KPEEDNDYHRSDEQVSICLE CNSSKLRGLKRKWIRCSAQA TVLHLKKFIKKLNLSSFNE LDILCNEEILGKDHTLKFFV VTRWRFKKAPLLHYPKM
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-242 of human PCGF3 (NP_006306.2).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of extracts of various cell lines, using PCGF3 antibody (STJ29595) at 1:1000 dilution. _Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. _Lysates/proteins: 25 Mu g per lane. _Blocking buffer: 3% nonfat dry milk in TBST. _Detection: ECL Enhanced Kit. _Exposure time: 30s.