

## Anti-TSSK3 antibody (1-268) (STJ29742)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Mouse/Rat

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 30kDa Observed Mw: 30kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">81629</a>
Gene Symbol	<a href="#">TSSK3</a>
UniProt ID	<a href="#">TSSK3_HUMAN</a>
Immunogen Region	1-268
Immunogen Sequence	MEDFLLSNGYQLGKTIGEGT YSKVKEAFSKKHQRKVAIKV IDKMGGPEEFIQRFLEPRELQ IVRTLDHKNIQVYEMLESA DGKICLVMELAEGGDVFDV LGGPLPESRAKALFROMVE AIRYCHGCGVAHRDLKCENA LLQGFNLKLTDFGFAKVLPK SHRELSQTFCGSTAYAAPEV LQGIPHDSKKGVDVWSMGVVL YVMLCASLPFDDTDIPKMLW QQKQGVSPFTHLSISADCQ
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-268 of human TSSK3 (NP_443073.1).

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Western blot analysis of various lysates using TSSK3 Rabbit polyclonal antibody (STJ29742) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.