

Anti-HDAC4 antibody (530-630) (STJ29788)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA: Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 119kDa Observed Mw: 140kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9759
Gene Symbol	HDAC4
UniProt ID	HDAC4_HUMAN
Immunogen Region	530-630
Immunogen Sequence	REHQALLDEPYLDRLPQGKE AHAQAGVQVKQEPIDESDEEE AEPPREVEPGQRQPSEQELL FRQQALLLEQQRIHQLRNYQ ASMEAAGIPVSFGGHRPLSR A
Specificity	A synthetic peptide corresponding to a sequence within amino acids 530-630 of human HDAC4 (NP_006028.2).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of various lysates, using HDAC4 Rabbit polyclonal antibody (STJ29788) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.