

Anti-hCG beta antibody [5012] (STJ400015)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	CLIA/FIA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	37 mM citrate, 125 mM phosphate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Antibody recognizes human chorionic gonadotropin free beta subunit. Does not recognize intact hCG, hCGa, LH, FSH, or TSH
-------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---