

## Anti-C-peptide antibody (4-19aa) [9103] (STJ400018)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FIA
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	37 mM citrate, 125 mM phosphate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage Instruction	Store at 2–8 °C, and do not freeze.

### TARGET INFORMATION

Immunogen Region	4-19aa
Specificity	Antibody recognizes human C-peptide. Human proinsulin (recombinant) 18 %

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---