

Anti-CA125 antibody [4601] (STJ400019)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/FIA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	50 mM citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Antibody recognizes human cancer antigen 125 (CA125). Does not recognize carbohydrate antigen 19-9 (CA19-9) , cancer antigen 15-3 (CA15-3) , CYFRA 21-1, carcinoembryonic antigen (CEA) , alpha fetoprotein (AFP) , or human epididymis protein 4 (HE4)
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---