

Anti-CA19-9 antibody [4701] (STJ400025)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FIA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	50 mM citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG3
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Antibody recognizes carbohydrate antigen 19-9 (CA19-9). Does not recognize cancer antigens CA125 or CA15-3. Does not recognize sialylated Lewis x (sLex) or Lewis y carbohydrate structures. Minor cross-reactivity against sialylated Lewis c (sLec) carb
-------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---