

Anti-Chlamydia trachomatis antibody [6701] (STJ400036)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FIA
Host / Source	Mouse
Reactivity	C.trachomatis

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Antibody recognizes Chlamydia trachomatis LPS. Recognizes C. pneumoniae, C. psittaci and S. aureus (unspecific binding to protein A). Cross-reactivity with the following micro-organisms below 5 %: A. lwoffii, A. baumannii, B. fragilis, C. albicans, C
-------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---