

## Anti-Chlamydia trachomatis antibody [6709] (STJ400038)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FIA
Host / Source	Mouse
Reactivity	C.trachomatis

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	0.9 % NaCl, 0.095 % NaN <sub>3</sub> as a preservative
Isotype	IgG2a
Storage Instruction	Store at 2–8 °C, and do not freeze.

### TARGET INFORMATION

Specificity	Antibody recognizes Chlamydia trachomatis LPS. Recognizes C. pneumoniae and C. psittaci. Does not recognize S. minnesota, B. quintana, B. henselae or S. negevensis.
-------------	--

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---