

**Anti-CK-MB antibody [7502] (STJ400040)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FIA                |
| Host / Source | Mouse              |
| Reactivity    | Human              |

**PRODUCT PROPERTIES**

|                     |                                                                                     |
|---------------------|-------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                          |
| Concentration       | 5.0 mg/ml (+/-10 %)                                                                 |
| Conjugation         | Unconjugated                                                                        |
| Formulation         | 37 mM citrate, 125 mM phosphate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative |
| Isotype             | IgG1                                                                                |
| Storage Instruction | Store at 2–8 °C, and do not freeze.                                                 |

**TARGET INFORMATION**

|             |                                                                         |
|-------------|-------------------------------------------------------------------------|
| Specificity | Antibody recognizes human creatine kinase MB isoenzyme. CK-BB isoenzyme |
|-------------|-------------------------------------------------------------------------|

**ADDITIONAL INFORMATION**

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|