

## Anti-CRP antibody [6404] (STJ400048)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/FIA/IT
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	> 6 mg/ml
Conjugation	Unconjugated
Formulation	0.9 % NaCl, 0.095 % NaN <sub>3</sub> as a preservative
Isotype	IgG1
Storage Instruction	Store at 2–8 °C, and do not freeze.

### TARGET INFORMATION

Specificity	Antibody recognizes human C-reactive protein (CRP). Does not recognize serum amyloid-P complex
-------------	--

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---