

## Anti-CYFRA21-1 antibody [1605] (STJ400060)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FIA                |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal  |
| Concentration       | 5.0 mg/ml (+/-10 %)   |
| Conjugation         | Unconjugated  |
| Formulation         | 50 mM citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN <sub>3</sub> as a preservative |
| Isotype             | IgG1  |
| Storage Instruction | Store at 2–8 °C, and do not freeze.   |

### TARGET INFORMATION

|             |   |
|-------------|---|
| Specificity | Antibody recognizes human CYFRA 21-1 fragment of cytokeratin 19. 100 % cytokeratin 19 |
|-------------|---|

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|