

Anti-hCG beta antibody [5009] (STJ400084)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	CLIA/FIA
Host / Source	Mouse
Reactivity	hCG beta

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	0.9 % NaCl, 0.095 % NaN ₃ as a preservative
Isotype	IgG1
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Antibody recognizes human chorionic gonadotropin. LH (KD = 1.2 nM). Does not recognize hCG α , FSH, or TSH
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---