

## Anti-hCoV OC43 antibody [3505] (STJ400092)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FIA
Host / Source	Mouse
Reactivity	hCoV OC43

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	50 mM citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN <sub>3</sub> as a preservative
Isotype	IgG1
Storage Instruction	Store at 2–8 °C, and do not freeze.

### TARGET INFORMATION

Specificity	Antibody recognizes human coronaviruses OC43 nucleoprotein. Does not recognize human coronaviruses 229E, NL63 or HKU1. Does not recognize recombinant SARS-CoV or SARS-CoV-2 nucleoprotein. Others not tested.
-------------	--

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---