

Anti-HRP2 antibody [3902] (STJ400102)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FIA
Host / Source	Mouse
Reactivity	HRP2

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	50 mM citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Antibody recognizes Plasmodium falciparum histidine-rich protein 2 (HRP2). Not determined
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---