

Anti-Influenza A antibody [7305] (STJ400118)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FIA
Host / Source	Mouse
Reactivity	Influenza A

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	37 mM citrate, 125mM phosphate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG2a
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Antibody recognizes Influenza A nucleoprotein (np). Recognizes Influenza A (H3N2) group 3C.2a1b and group 3C.3a, Influenza A (H1N1) pdm09 group 6B.1A1 and group 6B.1A5. Antibody also recognizes H5N3, H7N3, H9N2, and H5N1. NA Does not recognize Influe
-------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---