

Anti-Influenza B antibody [9906] (STJ400123)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | LF/FIA |
| Host / Source | Mouse |
| Reactivity | Influenza B |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Concentration | 5.0 mg/ml (+/-10 %) |
| Conjugation | Unconjugated |
| Formulation | 50 mM citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative |
| Isotype | IgG2b |
| Storage Instruction | Store at 2–8 °C, and do not freeze. |

TARGET INFORMATION

| | |
|-------------|--|
| Specificity | Antibody recognizes Influenza B nucleoprotein (np). Recognizes Influenza B/Yamagata clade 3, Influenza B/Victoria clade 1A, clade 1A (?2) , and clade 1A (?3). Does not recognize Influenza A (H3N2) group 3C.2a1b and group 3C.3a, Influenza A (H1N1) pdm |
|-------------|--|

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|