

Anti-LH antibody [5301] (STJ400126)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/CLIA/FIA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Antibody recognizes human luteinizing hormone (lutropin) , and its beta-subunit. LH a 13 %, LH β 170 %, FSH 5 %, hCG 138 %, TSH 0.03 %
-------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---