

Anti-LH antibody [5302] (STJ400127)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/CLIA/FIA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Antibody recognizes human luteinizing hormone (lutropin) , and its beta-subunit. LH a 10 %, LH β 127 %, FSH 3 %, hCG 4 %, TSH 0.02 %
-------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---