

## Anti-PG II antibody [8101] (STJ400157)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | LF/FIA/IT          |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                                   |
|---------------------|-------------------------------------------------------------------|
| Clonality           | Monoclonal                                                        |
| Concentration       | 5.0 mg/ml (+/-10 %)                                               |
| Conjugation         | Unconjugated                                                      |
| Formulation         | 50 mM citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative |
| Isotype             | IgG1                                                              |
| Storage Instruction | Store at 2–8 °C, and do not freeze.                               |

### TARGET INFORMATION

|             |                                                                                       |
|-------------|---------------------------------------------------------------------------------------|
| Specificity | Antibody recognizes human pepsinogen II. Human Pepsinogen I 2.2 %, Porcine Pepsinogen |
|-------------|---------------------------------------------------------------------------------------|

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|