

Anti-Progesterone antibody [1805] (STJ400165)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FIA
Host / Source	Mouse
Reactivity	Progesterone

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	37 mM citrate, 125 mM phosphate, pH 6.0, 0.9 % NaCl, 0.05 % Sulfobetaine, 0.095 % NaN3 as a preservative
Isotype	IgG2a
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Antibody recognizes progesterone. 11-alpha-hydroxyprogesterone 14 %, 17-alpha-hydroxyprogesterone 4 %, 21-hydroxyprogesterone 1 %, 17-alpha-hydroxypregnenolone 0 %
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---