

## Anti-TSH antibody [5407] (STJ400195)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | LF/CLIA/FIA        |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                                   |
|---------------------|-------------------------------------------------------------------|
| Clonality           | Monoclonal                                                        |
| Concentration       | 5.0 mg/ml (+/-10 %)                                               |
| Conjugation         | Unconjugated                                                      |
| Formulation         | 50 mM citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative |
| Isotype             | IgG1                                                              |
| Storage Instruction | Store at 2–8 °C, and do not freeze.                               |

### TARGET INFORMATION

|             |                                                                        |
|-------------|------------------------------------------------------------------------|
| Specificity | Antibody recognizes human thyrotropin and its free beta subunit. TSH a |
|-------------|------------------------------------------------------------------------|

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|