

## Anti-FABP3 antibody [2304] (STJ400204)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/FIA/IT
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	50 mM citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN <sub>3</sub> as a preservative
Isotype	IgG1
Storage Instruction	Store at 2–8 °C, and do not freeze.

### TARGET INFORMATION

Specificity	Antibody recognizes human fatty acid binding protein 3 (FABP3). Does not recognize FABP1 or FABP4.
-------------	--

### ADDITIONAL INFORMATION

Note	<b>STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.</b>
------	--