

Anti-cTnT antibody [3701] (STJ400209)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	Greater than 95% purity
Formulation	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage Instruction	Unspecified, storage at 2–8°C. Temperature analysis:-70 °C, 21 days OK-20 °C, 21 days OK +4 °C, 21 days OK +35 °C, 7 days OK +35 °C, 21 days OK +45 °C, 3 days OK +45 °C, 7 days OK. Shelf life stability testing at 2–8 °C in the product buffer.

TARGET INFORMATION

Gene ID	7139
Gene Symbol	TNNT2
UniProt ID	TNNT2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---