

Anti-ST2 antibody [10202] (STJ400211)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	Greater than 95% purity
Formulation	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage Instruction	Unspecified, storage at 2–8°C. Temperature analysis:-70 °C, 21 days OK-20 °C, 21 days OK +4 °C, 21 days OK +35 °C, 7 days OK +35 °C, 21 days Minor charge alterations +45 °C, 3 days OK +45 °C, 7 days Minor charge alterations. Shelf life stability test

TARGET INFORMATION

Gene ID	6822
Gene Symbol	SULT2A1
UniProt ID	ST2A1_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---