

## Anti-cTnT antibody (139-152) [3711] (STJ400222)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	Greater than 95% purity
Formulation	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage Instruction	Unspecified, storage at 2–8°C. Temperature analysis:-70 °C, 21 days Not Determined (N/D)-20 °C, 14 days OK +4 °C, 14 days OK +35 °C, 14 days OK +45 °C, 7 days OK. Shelf life stability testing at 2–8 °C in the product buffer.

### TARGET INFORMATION

Gene ID	<a href="#">7139</a>
Gene Symbol	<a href="#">TNNT2</a>
UniProt ID	<a href="#">TNNT2_HUMAN</a>
Immunogen Region	139-152

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---