

Anti-Galectin-3 antibody (66-104aa) [10305] (STJ400242)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	CLIA/FIA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	50 mM citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
Isotype	IgG1
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Immunogen Region	66-104aa
Specificity	Antibody recognizes human galectin-3. Does not recognize Galectin-1, Galectin-9, or Galectin-3-binding protein

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---