

## Anti-Vitamin 25 OH D2/D3 antibody [3851] (STJ400269)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Mouse
Reactivity	Vitamin 25

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	> 1.0 mg/ml
Conjugation	Unconjugated
Formulation	0.1 M PBS, pH 7.4, 0.9 % NaCl, 0.05 % NaN <sub>3</sub> as a preservative
Isotype	IgG1
Storage Instruction	Store at 2–8 °C, and do not freeze.

### TARGET INFORMATION

Specificity	Antibody recognizes natural 25 OH Vitamin D2 and 25 OH Vitamin D3. Vitamin D3
-------------	---

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---