

## Anti-Renin antibody [10851] (STJ400274)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA              |
| Host / Source | Mouse              |
| Reactivity    | Renin              |

### PRODUCT PROPERTIES

|                     |                                     |
|---------------------|-------------------------------------|
| Clonality           | Monoclonal                          |
| Concentration       | > 1.0 mg/ml                         |
| Conjugation         | Unconjugated                        |
| Formulation         | 0.1 M BBS, pH 8, 0.9 % NaCl         |
| Isotype             | IgG1                                |
| Storage Instruction | Store at 2–8 °C, and do not freeze. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">5972</a>       |
| Gene Symbol | <a href="#">REN</a>        |
| UniProt ID  | <a href="#">RENI_HUMAN</a> |
| Specificity | Human                      |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|