

Anti-Renin antibody (STJ400275)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Mouse
Reactivity	Renin

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	> 1.0 mg/ml
Conjugation	Unconjugated
Formulation	0.1 M PBS, pH 7.4, 0.9 % NaCl, 0.05 % NaN3 as a preservative
Isotype	IgG1
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Gene ID	5972
Gene Symbol	REN
UniProt ID	RENI_HUMAN
Specificity	Human

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---