

Anti-CEA antibody [5905] (STJ400281)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	CLIA/FIA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage Instruction	Store at 2-8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Antibody recognizes human carcinoembryonic antigen. Some cross-reactivity with human NCA (non-specific cross-reacting antigen) and NCA-2 (non-specific cross-reacting antigen-2)
-------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---