

Anti-cTnI antibody (190-195aa) [RC9707] (STJ400284)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FIA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	50 mM citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
Isotype	IgG1
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Immunogen Region	190-195aa
Specificity	Antibody recognizes human cardiac troponin I (cTnI). Skeletal troponin I (sTnI) 14 % Does not recognize cardiac troponin T (cTnT) , cardiac troponin C (cTnC) , or skeletal troponin T (sTnT)

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---