

## Anti-Estriol antibody [RC10601] (STJ400338)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FIA                |
| Host / Source | Mouse              |
| Reactivity    | E3                 |

### PRODUCT PROPERTIES

|                     |                                                                   |
|---------------------|-------------------------------------------------------------------|
| Clonality           | Monoclonal                                                        |
| Concentration       | 5.0 mg/ml (+/-10 %)                                               |
| Conjugation         | Unconjugated                                                      |
| Formulation         | 50 mM citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative |
| Isotype             | IgG1                                                              |
| Storage Instruction | Store at 2–8 °C, and do not freeze.                               |

### TARGET INFORMATION

|             |                                                                 |
|-------------|-----------------------------------------------------------------|
| Specificity | Antibody recognizes estriol. Does not recognize Estradiol (E2). |
|-------------|-----------------------------------------------------------------|

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|