

## Anti-ACTH antibody [R12801] (STJ400340)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	Liquid, may turn slightly opaque during storage
Isotype	IgG1

### TARGET INFORMATION

Specificity	Antibody recognizes human adrenocorticotrophic hormone
-------------	--

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---