

## Anti-ACTH antibody [R12802] (STJ400341)

### GENERAL INFORMATION

|              |                    |
|--------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ELISA              |
| Reactivity   | Human              |

### PRODUCT PROPERTIES

|               |                                                 |
|---------------|-------------------------------------------------|
| Clonality     | Monoclonal                                      |
| Concentration | 5.0 mg/ml (+/-10 %)                             |
| Conjugation   | Unconjugated                                    |
| Formulation   | Liquid, may turn slightly opaque during storage |
| Isotype       | IgG1                                            |

### TARGET INFORMATION

|             |                                                        |
|-------------|--------------------------------------------------------|
| Specificity | Antibody recognizes human adrenocorticotrophic hormone |
|-------------|--------------------------------------------------------|

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|